

Product list
2020

Product list of the PCR diagnostic kits
Compliance with the requirements of international standard EN ISO 13485:2016

Contents

Nucleic acids extraction	5
Blood-borne infections	7
Human immunodeficiency virus	7
Hepatitis viruses A, B, C, D, G.....	7
Quantification of human DNA in clinical specimens.....	9
Sexually transmitted infections (STIs).....	9
Vaginal dysbiosis	10
Mycoses	11
Human papillomaviruses	11
Human herpes viruses and TORCH-infections	12
Tick-borne infections	13
Nosocomial infections.....	14
Respiratory infections	14
Gastrointestinal infections.....	14
Tuberculosis	15
Single nucleotide polymorphisms (SNPs) and mutations in human genes.....	15

New kits

RealBest Sorbitus (variant 4x24)	6
RealBest Sorbitus (variant 1x96)	6
RealBest Bioflor	11
RealBest DNA Staphylococcus aureus/mecA/lukS-PV (set 1)	14
RealBest DNA Staphylococcus aureus/mecA/lukS-PV (set 2)	14
RealBest RNA SARS-CoV-2	14

The RealBest diagnostic kits are intended for the detection of infectious agents and determination of genetic polymorphisms by real-time PCR. In-house developed technologies turn the procedure of real-time PCR analysis from a complicated process requiring extensive specialized knowledge and practical experience into a standardized, simple, and readily reproducible procedure. AO Vector-Best offers user-friendly solutions for diagnostics such as lyophilized ready-to-use mixtures with a long shelf-life and transportation temperature of up to 25 °C.

RealBest kits:

- Complex approach to PCR diagnostics encompassing sample preparation, PCR and evaluation of results
- Simple nucleic acids extraction using magnetic particles
- Lyophilized Ready Master Mix for the detection of DNA/RNA in real-time
- Compatible PCR and RT-PCR protocols for the detection of a wide range of pathogens
- Noncompetitive Internal Control
- Storage of all components at the same temperature (2–8)°C for up to 12 months; no freezing required
- Transportation at the temperature up to 25 °C for no more than 10 days is acceptable

Different formats:

Str-format: Ready Master Mix is lyophilized in 0.2 ml white tubes of the 96-tube plate. The only step required for the reaction procedure is the addition of extracted DNA/RNA samples.

Fla-format: Master Mix for 10 reactions is lyophilized in 2.0 ml tubes. After dilution with Recovery Solution, dispense Master Mix to the optimized tubes, add the extracted DNA/RNA, and run the test.

Uni-format: Ready Master Mix is lyophilized in clear 0.2 ml tubes. Add the extracted DNA/RNA into the tubes and run the test. One tube – one reaction.

Real-time PCR Systems and formats:





Str, Uni and **Fla** formats can be used with: CFX96™ (Bio-Rad, USA), DT*prime* (DT-96), DT*lite* (DNA-Technology, Russia).

Fla, Uni can be used with: Rotor-Gene® Q (Qiagen, Germany), Rotor-Gene® 6000/3000 (Corbett Research, Australia).

Number of simultaneously detected pathogens	Detection channels
1	FAM, ROX
2	FAM, ROX, HEX (JOE)
3 and more	FAM, ROX, HEX (JOE), Cy5

Other real-time PCR cyclers with similar technical characteristics can be used after validation.

Nucleic acids extraction

Cat. No.	Kit name	Number of tests	Clinical specimens	Comments
Manual nucleic acids extraction				
C-8891	RealBest DeltaMag HBV/HCV/HIV (variant 1-8)	48 (6x8)	Blood serum/plasma Sample volume – 1 ml	Extraction of HBV DNA HCV RNA HIV RNA
C-8893	RealBest DeltaMag HBV/HCV/HIV (variant 0.25-8)	48 (6x8)	Blood serum/plasma Sample volume – 250 µl	Extraction of HBV DNA HCV RNA HIV RNA
C-8895	RealBest extraction 1000	48 (6x8)	Blood serum/plasma Sample volume – 1 ml	Extraction of DNA and RNA
C-8896 	RealBest extraction 100	48 (6x8)	Whole blood, blood serum/plasma, leukocyte blood fraction, biopsy material, cerebrospinal fluid, urine, feces, epithelial cell specimens, tick suspensions, water samples Sample volume – 100 µl	Extraction of DNA and RNA
C-8889 	RealBest DNA-extraction 3 (variant 6x16)	96 (6x16)	Blood serum/plasma, urine, epithelial cell specimens Sample volume – 100 µl	Extraction of DNA
C-8897 	RealBest DNA-extraction 2	96 (4x24)	Blood serum/plasma, biopsy material, cerebrospinal fluid, urine, epithelial cell specimens Sample volume – 100 µl	Extraction of DNA
C-8899 	RealBest DNA-express	100	Epithelial cell specimens from the cervical canal, urethra, vagina, posterior laryngeal wall, conjunctiva; semen, saliva, urine	Extraction of DNA in 15 minutes. Lysis Reagent – 100 tubes, Transport solution – 100 tubes
C-8845	RealBest-Genetics DNA- express	50	Whole blood, buccal epithelium	Extraction of DNA in 20 minutes
Automated nucleic acids extraction				
C-8849	RealBest DeltaMag HBV/HCV/HIV (variant 1-16)	48 (3x16)	Blood serum/plasma Sample volume – 1 ml	Extraction of HBV DNA HCV RNA HIV RNA (for robotic platforms TECAN Freedom EVO)
C-8850	RealBest DeltaMag HBV/HCV/HIV (variant 0.25-16)	48 (3x16)	Blood serum/plasma Sample volume – 250 µl	Extraction of HBV DNA HCV RNA HIV RNA (for robotic platforms TECAN Freedom EVO)

Cat. No.	Kit name	Number of tests	Clinical specimens	Comments
C-8880	RealBest DeltaMag HBV/HCV/HIV (variant 0.25-48)	96 (2x48)	Blood serum/plasma Sample volume – 250 µl	Extraction of HBV DNA HCV RNA HIV RNA (for KingFisher Flex systems)
C-8847	RealBest Sorbitus variant 4x24	96 (4x24)	Nasal, nasopharyngeal, and oropharyngeal swabs, sputum, bronchoalveolar lavage fluid Recommended for respiratory infections Sample volume – 100 µl	Extraction of DNA and RNA (for any open-type automated NA extraction systems)
C-8848	variant 1x96	96 (1x96)		
C-8883 CE	RealBest UniMag	96 (4x24)	Whole blood, blood serum/plasma, leukocyte blood fraction, biopsy material, cerebrospinal fluid, urine, feces, epithelial cell specimens, tick suspensions, water samples Sample volume – 100 µl	Extraction of DNA and RNA (for KingFisher Flex systems)
C-8878	RealBest DNA-extraction 3 (variant 2x48)	96 (2x48)	Blood serum/plasma, urine, epithelial cell specimens Sample volume – 100 µl	Extraction of DNA (for KingFisher Flex systems and robotic platforms TECAN Freedom EVO)
Solutions and components				
C-8814	RealBest Hemolytic	50	Whole blood	Solution for Hemolysis 2 vials, 14 ml each
C-8879	Lysis Reagent (A)	192	Epithelial cell specimens from the cervical canal, urethra, vagina, posterior laryngeal wall, conjunctiva; semen, saliva, urine	Extraction of DNA 1 vial, 140 ml
C-8887	Lysis Reagent (M)	200	Epithelial cell specimens from the cervical canal, urethra, vagina, posterior laryngeal wall, conjunctiva; semen, saliva, urine	Extraction of DNA 200 tubes, 500 µl each
C-8894 CE	Transport solution (1)	100	Biopsy material, epithelial cell specimens from the cervical canal, urethra, vagina, posterior laryngeal wall, conjunctiva	100 tubes, 300 µl each
C-8885 CE	Transport solution (2)	200	Biopsy material, epithelial cell specimens from the cervical canal, urethra, vagina, posterior laryngeal wall, conjunctiva	200 tubes, 1000 µl each
C-8881 CE	Internal Control sample	192	Control of the efficiency of nucleic acids extraction and the absence of PCR inhibitors in specimens	Kit contents: – Internal Control sample, lyophilized – 6 vials; – Recovery Solution for Control samples – 2 vials, 4 ml each; – Specimen Diluent - 4 vials, 15 ml each

Blood-borne infections

<i>Cat. No.</i>	<i>Kit name</i>	<i>Format</i>	<i>Number of tests</i>	<i>Sensitivity</i>	<i>Comments</i>
Differential detection of human immunodeficiency virus types 1 and 2 RNA, hepatitis B virus DNA, hepatitis C virus RNA					
D-0592	RealBest HBV/HCV/HIV PCR	Str	48	HBV - 10 IU/ml HCV - 15 IU/ml HIV-1 - 30 IU/ml HIV-2 - 50 IU/ml	Qualitative detection


Human immunodeficiency virus

<i>Cat. No.</i>	<i>Kit name</i>	<i>Format</i>	<i>Number of tests</i>	<i>Sensitivity</i>	<i>Comments</i>
Proviral DNA of human immunodeficiency virus (HIV-1)					
D-0186	RealBest DNA HIV (WB)	Uni*	48	25 copies in the sample volume	Qualitative detection Clinical specimen: whole blood
D-0187	RealBest DNA HIV (DS)	Uni*	48	50 copies in the sample volume	Qualitative detection Clinical specimen: dry blood stains

Human immunodeficiency virus RNA


D-0198	RealBest RNA HIV Form 1	Uni*	48	20 IU/ml	Qualitative detection
D-0197	Form 2	Uni	48		
D-0199	RealBest RNA HIV quantitative	Uni*	48	20 IU/ml	Quantification
D-0195	RealBest RNA HIV Set 1	Str	96	20 IU/ml	Quantification
D-0196	Set 2	Str	96		Qualitative detection

Hepatitis viruses A, B, C, D, G

<i>Cat. No.</i>	<i>Kit name</i>	<i>Format</i>	<i>Number of tests</i>	<i>Sensitivity</i>	<i>Comments</i>
Hepatitis A virus RNA					
 D-0398	RealBest RNA HAV Set 1	Str	48	50 copies in the sample	Qualitative detection
D-0397	Set 2	Fla	50		

Hepatitis B virus DNA

D-0598	RealBest DNA HBV Form 1	Uni*	48	5 IU/ml	Qualitative detection
D-0597	Form 2	Uni	48		
D-0599	RealBest DNA HBV quantitative	Uni*	48	5 IU/ml	Quantification

<i>Cat. No.</i>	<i>Kit name</i>	<i>Format</i>	<i>Number of tests</i>	<i>Sensitivity</i>	<i>Comments</i>
	RealBest HBV PCR				
D-0595	Variant 1/quantitative	Str	96	5 IU/ml	Quantification
D-0596	Variant 2	Str	96		Qualitative detection
Hepatitis C virus RNA					
	RealBest RNA HCV				
D-0798	Form 1	Uni*	48	15 IU/ml	Qualitative detection
D-0790	Form 2	Uni	48		
D-0799	RealBest RNA HCV quantitative	Uni*	48	15 IU/ml	Quantification
	RealBest RNA HCV				
D-0794	Set 1	Str	96	15 IU/ml	Quantification
D-0795	Set 2	Str	96		Qualitative detection
Genotyping of hepatitis C virus					
D-0797	RealBest RNA HCV-genotype	Str	48	15 IU/ml	Quantification and genotyping
				400 IU/ml	genotyping
D-0793	RealBest RNA HCV-1/2/3	Uni	48	400 IU/ml	Genotyping
Hepatitis D virus RNA					
	RealBest RNA HDV				
D-0998	Set 1	Str	48	10 copies/ml	Qualitative detection
D-0997	Set 2	Fla	50		
Hepatitis G virus RNA					
	RealBest RNA HGV				
D-1298	Set 1	Str	48	50 copies/ml	Qualitative detection
D-1299	Set 2	Fla	50		

Kits **Uni*** contain reagents for nucleic acids extraction

Quantification of human DNA in clinical specimens

<i>Cat. No.</i>	<i>Kit name</i>	<i>Format</i>	<i>Number of tests</i>
CE	RealBest Sample Validation		
D-8888	Set 1	Str	96
D-8884	Set 2	Fla	100
Sexually transmitted infections (STIs)			
<i>Cat. No.</i>	<i>Kit name</i>	<i>Format</i>	<i>Number of tests</i>
CE0123	RealBest DNA Chlamydia trachomatis		
D-1998	Set 1	Str	96
D-1996	Set 2	Fla	100
CE	RealBest DNA Trichomonas vaginalis		
D-2098	Set 1	Str	96
D-2096	Set 2	Fla	100
CE	RealBest DNA Mycoplasma genitalium		
D-4396	Set 1	Str	96
D-4394	Set 2	Fla	100
CE	RealBest DNA Neisseria gonorrhoeae		
D-4498	Set 1	Str	96
D-4496	Set 2	Fla	100
CE	RealBest DNA Neisseria gonorrhoeae test 2		
D-4494	Set 1	Str	96
D-4495	Set 2	Fla	100
CE	RealBest DNA Treponema pallidum		
D-1898	Set 1	Str	96
D-1896	Set 2	Fla	100
Multiplex kits for STIs detection			
	RealBest DNA Chlamydia trachomatis/Ureaplasma species		
D-0455	Set 1	Str	96
D-0456	Set 2	Fla	100
CE0123	RealBest DNA Chlamydia trachomatis/Ureaplasma urealyticum		
D-0492	Set 1	Str	96
D-0491	Set 2	Fla	100
CE0123	RealBest DNA Chlamydia trachomatis/Mycoplasma genitalium		
D-0490	Set 1	Str	96
D-0487	Set 2	Fla	100
CE0123	RealBest DNA Chlamydia trachomatis/Trichomonas vaginalis		
D-0498	Set 1	Str	96
D-0497	Set 2	Fla	100
CE0123	RealBest DNA Chlamydia trachomatis/Neisseria gonorrhoeae		
D-0457	Set 1	Str	96
D-0458	Set 2	Fla	100
CE	RealBest DNA Mycoplasma hominis/Mycoplasma genitalium		
D-0494	Set 1	Str	96
D-0493	Set 2	Fla	100

Cat. No.	Kit name	Format	Number of tests
CE D-0496	RealBest DNA Trichomonas vaginalis/Neisseria gonorrhoeae Set 1	Str	96
D-0495	Set 2	Fla	100
CE D-0477	RealBest DNA Trichomonas vaginalis/Gardnerella vaginalis Set 1	Str	96
D-0478	Set 2	Fla	100
CE D-0444	RealBest DNA Candida albicans/Gardnerella vaginalis Set 1	Str	96
D-0445	Set 2	Fla	100
D-0488	RealBest PCR-12 STD*	Str	96

Vaginal dysbiosis

Cat. No.	Kit name	Format	Number of tests
CE D-2292	RealBest DNA Ureaplasma species Set 1	Str	96
D-2293	Set 2	Fla	100
CE D-2298	RealBest DNA Ureaplasma urealyticum Set 1	Str	96
D-2296	Set 2	Fla	100
CE D-4398	RealBest DNA Mycoplasma hominis Set 1	Str	96
D-4393	Set 2	Fla	100
CE D-4598	RealBest DNA Gardnerella vaginalis Set 1	Str	96
D-4596	Set 2	Fla	100
D-4205	RealBest Lactonorm** Set 1	Str	96
D-4206	Set 2	Fla	100
D-4228	RealBest DNA Streptococcus agalactiae Set 1	Str	96
D-4229	Set 2	Fla	100

Multiplex kits for dysbiosis detection

CE D-2294	RealBest DNA Ureaplasma urealyticum/Ureaplasma parvum Set 1	Str	96
D-2295	Set 2	Fla	100
D-4207	RealBest DNA Gardnerella vaginalis/Atopobium vaginae Set 1	Str	96
D-4208	Set 2	Fla	100
D-4212	RealBest DNA Prevotella species/Leptotrichia amnionii group Set 1	Str	96
D-4213	Set 2	Fla	100

* *Chlamydia trachomatis*, *Ureaplasma species*, *Mycoplasma hominis*, *Mycoplasma genitalium*, *Trichomonas vaginalis*, *Neisseria gonorrhoeae*, *Candida albicans*, *Gardnerella vaginalis*, cytomegalovirus, herpes simplex virus types 1 and 2, human papillomavirus types 16 and 18 (differential detection)

** Determination of the percentage of *Lactobacillus* sp. in the total number of bacteria

Cat. No.	Kit name	Format	Number of tests
CE	RealBest DNA Mobiluncus mulieris/Mobiluncus curtisii		
D-4201	Set 1	Str	96
D-4202	Set 2	Fla	100
CE	RealBest DNA Saccharimonas aalborgensis (TM7)/BVAB2		
D-4214	Set 1	Str	96
D-4215	Set 2	Fla	100
D-4225	RealBest Bioflor*	Str	96

Mycoses

Cat. No.	Kit name	Format	Number of tests
CE	RealBest DNA Candida albicans		
D-4698	Set 1	Str	96
D-4696	Set 2	Fla	100

Multiplex kits for mycoses detection

CE	RealBest DNA Candida albicans/Fungi		
D-0448	Set 1	Str	96
D-0449	Set 2	Fla	100
CE	RealBest DNA Candida parapsilosis/Candida tropicalis		
D-0440	Set 1	Str	96
D-0441	Set 2	Fla	100
CE	RealBest DNA Candida famata/Candida guilliermondii		
D-0442	Set 1	Str	96
D-0443	Set 2	Fla	100
CE	RealBest DNA Candida krusei/Candida glabrata		
D-0446	Set 1	Str	96
D-0447	Set 2	Fla	100

Human papillomaviruses

Cat. No.	Kit name	Format	Number of tests
	RealBest DNA HPV 44		
D-8447	Set 1	Str	96
D-8443	Set 2	Fla	100
CE	RealBest DNA HPV 66		
D-8448	Set 1	Str	96
D-8440	Set 2	Fla	100
CE	RealBest DNA HPV 68		
D-8488	Set 1	Str	96
D-8489	Set 2	Fla	100

* Determination of the percentage of *Lactobacillus sp.* and bacteria causing vaginosis and vaginitis (*Gardnerella vaginalis*, *Atopobium vaginae*, *Prevotella species*, *Leptotrichia amnionii group*, *Streptococcus species*, *Enterococcus species*, and *Staphylococcus species*) in the total number of bacteria

<i>Cat. No.</i>	<i>Kit name</i>	<i>Format</i>	<i>Number of tests</i>
Multiplex kits for HPV detection			
CE D-8475	RealBest DNA HPV 6/11 Set 1	Str	96
D-8481	Set 2	Fla	100
CE D-8473	RealBest DNA HPV 16/18 Set 1	Str	96
D-8474	Set 2	Fla	100
D-8459	RealBest DNA HPV 26/51 Set 1	Str	96
D-8457	Set 2	Fla	100
CE D-8471	RealBest DNA HPV 31/33 Set 1	Str	96
D-8480	Set 2	Fla	100
CE D-8469	RealBest DNA HPV 35/45 Set 1	Str	96
D-8477	Set 2	Fla	100
CE D-8446	RealBest DNA HPV 6/11/44 Set 1	Str	96
CE D-8449	RealBest DNA HPV 26/53/66 Set 1	Str	96
CE D-8451	RealBest DNA HPV 68/73/82 Set 1	Str	96
CE D-8444	RealBest DNA HPV HR screen (detection of human papillomavirus types 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68 DNA) Set 1	Str	96
D-8445	Set 2	Fla	100
CE D-8479	RealBest DNA HPV HR genotype (differential detection of human papillomavirus types 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59 DNA) Set 1	Str	96
D-8482	Set 2	Fla	99
CE D-8478	RealBest DNA HPV HR genotype quantitative (differential detection and quantification of human papillomavirus types 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59 DNA) Set 1	Str	96
D-8487	Set 2	Fla	99

Human herpes viruses and TORCH-infections

<i>Cat. No.</i>	<i>Kit name</i>	<i>Format</i>	<i>Number of tests</i>
D-1798	RealBest DNA Toxoplasma gondii Set 1	Str	48
D-1796	Set 2	Fla	50
CE₀₁₂₃ D-2598	RealBest RNA Rubella Set 1	Str	48
D-2596	Set 2	Fla	50
CE D-2801	RealBest DNA Parvovirus B19 Set 1	Str	96
CE D-2193	RealBest DNA HSV-1,2 Set 1	Str	96
D-2194	Set 2	Fla	100

<i>Cat. No.</i>	<i>Kit name</i>	<i>Format</i>	<i>Number of tests</i>
CE	RealBest DNA VZV		
D-2185	Set 1	Str	48
D-2187	Set 2	Fla	50
CE	RealBest DNA EBV (quantification is possible)		
D-2198	Set 1	Str	96
D-2196	Set 2	Fla	100
CE0123	RealBest DNA CMV (quantification is possible)		
D-1598	Set 1	Str	96
D-1596	Set 2	Fla	100
CE	RealBest DNA HHV-6		
D-2150	Set 1	Str	48
D-2151	Set 2	Fla	50
	RealBest DNA HHV-8		
D-2148	Set 1	Str	48
D-2149	Set 2	Fla	50
Multiplex kits for human herpes viruses detection			
CE	RealBest DNA HSV-1/HSV-2		
D-2195	Set 1	Str	96
D-2197	Set 2	Fla	100
	RealBest DNA CMV/HSV-1,2		
D-0489	Set 1	Str	96
D-0486	Set 2	Fla	100
Tick-borne infections			
<i>Cat. No.</i>	<i>Kit name</i>	<i>Format</i>	<i>Number of tests</i>
CE	RealBest RNA TBEV		
D-1198	Set 1	Str	48
D-1199	Set 2	Fla	50
CE	RealBest DNA Borrelia burgdorferi s.l.		
D-1498	Set 1	Str	48
D-1499	Set 2	Fla	50
CE	RealBest DNA Borrelia miyamotoi		
D-1495	Set 1	Str	48
D-1496	Set 2	Fla	50
	RealBest DNA Babesia species		
D-5389	Set 1	Str	48
D-5390	Set 2	Fla	50
	RealBest DNA Rickettsia species		
D-5391	Set 1	Str	48
D-5392	Set 2	Fla	50
Multiplex tick-borne infections detection			
CE	RealBest DNA Borrelia burgdorferi s.l./RNA TBEV		
D-5396	Set 1	Str	48
	RealBest DNA Rickettsia sibirica/Rickettsia heilongjiangensis		
D-5393	Set 1	Str	48
D-5394	Set 2	Fla	50
CE	RealBest DNA Anaplasma phagocytophilum/Ehrlichia muris, Ehrlichia chaffeensis		
D-5398	Set 1	Str	48
D-5399	Set 2	Fla	50

Nosocomial infections

<i>Cat. No.</i>	<i>Kit name</i>	<i>Format</i>	<i>Number of tests</i>
	RealBest DNA Enterococcus faecalis/Enterococcus faecium		
D-4222	Set 1	Str	96
D-4226	Set 2	Fla	100
	RealBest DNA Acinetobacter baumannii/Stenotrophomonas maltophilia		
D-5601	Set 1	Str	96
D-5602	Set 2	Fla	100
	RealBest DNA Staphylococcus aureus/mecA/lukS-PV		
D-5603	Set 1	Str	96
D-5604	Set 2	Fla	100
	RealBest DNA Klebsiella pneumoniae/Pseudomonas aeruginosa		
D-5605	Set 1	Str	96
D-5606	Set 2	Fla	100

Respiratory infections

<i>Cat. No.</i>	<i>Kit name</i>	<i>Format</i>	<i>Number of tests</i>
D-5580 CE	RealBest RNA SARS-CoV-2	Str	96
	RealBest DNA Bordetella species/Bordetella pertussis/Bordetella bronchiseptica		
D-5586	Set 1	Str	48
D-5587 CE	Set 2	Fla	50
	RealBest DNA Streptococcus pneumoniae		
D-5590	Set 1	Str	48
D-5591 CE	Set 2	Fla	50
	RealBest DNA Haemophilus influenzae		
D-5592	Set 1	Str	48
D-5593 CE	Set 2	Fla	50
	RealBest DNA Chlamydophila pneumoniae		
D-5594 CE ₀₁₂₃	Set 1	Str	48
D-5595	Set 2	Fla	50
	RealBest DNA Mycoplasma pneumoniae		
D-5596 CE	Set 1	Str	48
D-5597	Set 2	Fla	50
	RealBest DNA Legionella pneumophila		
D-5598	Set 1	Str	48
D-5599	Set 2	Fla	50

Gastrointestinal infections

<i>Cat. No.</i>	<i>Kit name</i>	<i>Format</i>	<i>Number of tests</i>
	RealBest DNA Helicobacter pylori		
D-3798 CE	Set 1	Str	48
D-3796	Set 2	Fla	50

Tuberculosis

Cat. No.	Kit name	Format	Number of tests
CE D-2398	RealBest DNA MBTC Form 1	Uni*	48
D-2399	Form 2	Str	96

Single nucleotide polymorphisms (SNPs) and mutations in human genes

Cat. No.	Kit name	Format	Number of tests	Gene: polymorphism
Thrombophilia and hyperhomocysteinemia				
D-3801	RealBest-Genetics Hemostasis (12)	Str	48	F2:20210G/A; F5:1691G/A; F7:10976G/A; F13A1: c.103G/T; FGB:-455G/A; ITGA2:807C/T; ITGB3:1565T/C; PAI-1:- 6755G/4G; MTR:2756 A/G; MTRR:66 A/G; MTHFR:677C/T; MTHFR:1298A/C
D-3802 CE	RealBest-Genetics Hemostasis (F2/F5)	Str	48	F2:20210G/A; F5:1691G/A
D-3803 CE	RealBest-Genetics Hemostasis (MTR/MTRR/MTHFR)	Str	48	MTR:2756A/G; MTRR:66A/G; MTHFR:677C/T; MTHFR:1298A/C
D-3831 CE	RealBest-Genetics Hemostasis FGB/F13A1	Str	48	FGB: -455G/A; F13A1: c.103G/T
D-3832 CE	RealBest-Genetics Hemostasis ITGA2/F7	Str	48	ITGA2:807C/T; F7:10976G/A
D-3833 CE	RealBest-Genetics Hemostasis PAI-1/ITGB3	Str	48	PAI-1:-6755G/4G; ITGB3:1565T/C
Metabolic disorder				
D-3804 CE	RealBest-Genetics MCM6	Str	48	MCM6: -13910C/T
Cardiovascular diseases				
D-3805	RealBest-Genetics NOS3 T(-786)C/G894T	Str	48	NOS3: T(-786)C; G894T
D-3806	RealBest-Genetics ACE Del287/AGT T704C	Str	48	ACE: Del287; AGT: T704C
D-3817	RealBest-Genetics CYP11B2 C(-344)T/ADD1 G1378T	Str	48	CYP11B2: C(-344)T; ADD1: G1378T
D-3818	RealBest-Genetics AGT C521T/AGTR1 A1166C	Str	48	AGT: C521T; AGTR1: A1166C
D-3819	RealBest-Genetics GNB3 C825T/AGTR2 G1675A	Str	48	GNB3: C825T; AGTR2: G1675A
Predicting the success of HCV therapy				
D-3811 CE	RealBest-Genetics Interleukin 28B	Str	48	IL28B: rs12979860 C/T; rs8099917 T/G

Kit **Uni*** contains reagents for nucleic acids extraction

Cat. No.	Kit name	Format	Number of tests	Gene: polymorphism
Reproductive function				
D-3814 CE	RealBest-Genetics AZF microdeletions	Str	48	AZFa locus: sY-84, sY-86 AZFb locus: sY-127, sY-134 AZFc locus: sY-254, sY-255
Hereditary breast and ovarian cancer				
D-3807	RealBest-Genetics BRCA1 185delAG/3875del4	Str	48	BRCA1: 185delAG; 3875del4
D-3808	RealBest-Genetics BRCA1 3819del5/T300G	Str	48	BRCA1: 3819del5; T300G
D-3809	RealBest-Genetics BRCA1 2080delA (insA)/BRCA2 6174delT	Str	48	BRCA1: 2080delA (insA); BRCA2: 6174delT
D-3812	RealBest-Genetics BRCA1 4153delA/5382insC	Str	48	BRCA1: 4153delA; 5382insC
D-3810	RealBest-Genetics NBS1	Str	48	NBS1: 657del5
D-3816	RealBest-Genetics CHEK2	Str	48	CHEK2: 1100delC; IVS2+1G>A
Hereditary HFE-associated hemochromatosis (hemochromatosis type I)				
D-3822 CE	RealBest-Genetics Hemochromatosis (HFE 187 C/G)	Str	48	HFE: 187 C/G (His63Asp)
D-3823 CE	RealBest-Genetics Hemochromatosis (HFE 193 A/T)	Str	48	HFE: 193 A/T (Ser65Cys)
D-3824 CE	RealBest-Genetics Hemochromatosis (HFE 845 G/A)	Str	48	HFE: 845 G/A (Cys282Tyr)
Autoimmune disease				
D-3836	RealBest-Genetics HLA-B*27	Str	48	HLA-B*27
Pharmacogenetics of warfarin				
D-3827 CE	RealBest-Genetics Warfarin (CYP2C9*2/CYP2C9*3)	Str	48	CYP2C9: 430 C/T; 1075 A/C
D-3828 CE	RealBest-Genetics Warfarin (VKORC1/CYP4F2*3)	Str	48	VKORC1: 1173 C/T; CYP4F2: 1297 C/T
D-3829 CE	RealBest-Genetics Warfarin (GGCX)	Str	48	GGCX: rs11676382 C/G
D-3830 CE	RealBest-Genetics Warfarin (CYP2C9*5/CYP2C9*6)	Str	48	CYP2C9: 1080 C/G; 818 delA